ALE-HOP

User manual | Manual de usuario | Manual do usuário
Mode d'emploi | Benutzerhandbuch | Manuale d'uso
Priručnik s uputama | للتحويل بين الأوضاع

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PRODUCT SPECIFICATIONS

VOLTAGE: 5V/1A

POWER: 2W

BATTERY CAPACITY: 350MAH

HOW TO USE

1. Switch machine:

To hold product by hand, it will be starting up if short press the power button. It will be shut down if short press the power button again.

2. Massage:

Put the massage head on body massage parts.

3. Notice to charging:

- a. Plug the charging cable into the charging plug of the product.
- b. During charging, the power indicator light flashes blue; after fully charged, the battery indicator is off.
- c. The product does not support charging under starting up state. Do not use the product during charging.
- d. It is a best charging time for 2-3 hours.

CLEANING & MAINTENANCE

1. The product is electronic product. Do not rinse the product with cleaning solution or water.
2. Do not clean the product with any steel mesh, copper brush,

hard objects, gasoline and kerosene and others.

3. Please keep it in storage box, and place in a dry cool place when you do not use it.

NOTICES

To ensure you are correct use of the product, please read the following notices:

1. Do not use the product at high temperature.
2. Please do not use the product, if a contact part of massage head is redness and inflammation.
3. Please do not use the product if you have haemorrhagic disorders.
4. This product is designed for home electronic health care, only use for personal or non-professional.
5. Do not charge in an environment of excessive humidity.
6. Do not allow children to use the product directly.

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CLAVE DENIA S.A.